

Japanese Standard Medical Gas Outlets



We are a Chinese manufacturer with a factory specializing in medical equipment. One of our most popular products is the DFTT Japanese Standard Medical Gas Outlets. The reason why Japanese Standard Gas Outlets are so popular is that we have high cost performance, superior workmanship, and ample stock that delivery is never delayed. In the external procurement, we will only buy high-quality raw materials, never shoddy. Therefore, we hope to be your partner and I hope you can give me this opportunity.

Product Description

DFTT® **Japanese Standard Medical Gas Outlets** is similar to other Gas terminals, are made of copper and of excellent quality, which can meet all the special needs of our customers with complete functions. Japanese Standard Gas Outlets adopt cutting-edge technology, all-metal material, not easy to damage, good sealing, not easy to water, is your best choice.

Product Parameter

Product name	Shape	material	Size	color
Japanese Standard Gas Outlets	irregular shape	Copper	7*6*6CM	Golden

Japanese Standard Gas Outlets Feature And Application





Japanese Standard Gas Outlets are more convenient than other types of products. With convenient connection of pipelines, different types of Gas can be distinguished automatically by different locations of the pin holes, thus preventing accidental insertion and contact, which greatly ensures safety. From the characteristics of the point of view, the thread after careful grinding, in strict accordance with the national standards to create, simple rotation, strong and stable. Japanese Standard Gas Outlets are often used at the bedside of hospital wards, which fit perfectly with the wards.

Japanese Standard Gas Outlets Details





DFTT® **Japanese Standard Medical Gas Outlets** built-in access valve, color optional, can avoid disassembly panel with air maintenance, easy installation and maintenance, single construction. Japanese Standard Gas Outlets are tough and durable in whole. Japanese Standard Gas Outlets can be pressed with one hand, which is the biggest selling point over other products. Japanese Standard Gas Outlets have passed the quality system standards, as well as the factory air tightness test and 20000 plug and pull fatigue tests.